



ResMed

Actual Size: 140 mm (Width)
x 175 mm (Height)

Ruler is for scaling purposes.

How to use the fitting template:



1. Please ensure the template is printed to scale
2. Cut out along dotted line
3. Refer to diagram on template to determine appropriate cushion size

AirFit™ F20

Full face mask
Full face mask **for Her**


© 2017 ResMed Ltd. All Rights Reserved.
Specifications may change without notice.
Product availability may vary across regions.
1031823/2 2017-03

CE
Rx Only 0123



AirFit™ F20

Full face mask
Full face mask **for Her**



© 2017 ResMed Ltd. All Rights Reserved. Specifications may change without notice.
Product availability may vary across regions.

Rx Only **CE**
0123



ResMed